Sonosite SII Ultrasound System

User Guide Errata

Manufacturer	EC Authorized Representative	Australia Sponsor
FUJIFILM SonoSite, Inc.	FUJIFILM SonoSite B.V.	FUJIFILM SonoSite Australasia Pty Ltd
21919 30th Drive SE	Joop Geesinkweg 140	114 Old Pittwater Road
Bothell, WA 98021 USA	1114 AB Amsterdam,	BROOKVALE, NSW, 2100
T: 1-888-482-9449 or 1-425-951-1200	The Netherlands	Australia
F: 1-425-951-1201		

Caution

United States federal law restricts this device to sale by or on the order of a physician.

Sonosite SII is a trademark and registered trademark of FUJIFILM SonoSite, Inc. in various jurisdictions. SONOSITE and the SONOSITE logo are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions.

FUJIFILM is a registered trademark of FUJIFILM Corporation in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P23649-01

Publication date: August 2021

Copyright © 2021 FUJIFILM SonoSite, Inc. All rights reserved.

Sonosite SII Ultrasound System User Guide Errata

Introduction	1
Document conventions	1
Getting help.....	2
Getting started	3
System controls.....	3
Intended uses.....	4
Imaging	4
2D imaging.....	4
ECG.....	4
Measurements and calculations	4
Cardiac calculations	4
Cleaning and disinfecting	5
Cleaning and disinfecting the ECG cable and slave cable	5
Safety	5
Electrical safety classification.....	5
Electrical safety.....	6
Compatible accessories and peripherals	6
Labeling symbols	7

Introduction

The content in this document corrects the Sonosite SII Ultrasound System user guide in the following chapters. The revisions will be made in the next update.

- ▶ Getting Started
- ▶ Imaging
- ▶ Measurements and Calculations
- ▶ Cleaning and Disinfecting
- ▶ Safety

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ♦.

For a description of labeling symbols that appear on the product, see “Labeling Symbols” in the user guide.

Getting help

FUJIFILM SonoSite Technical Support

United States or Canada	+1 (877) 657-8118
Europe and Middle East	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia and Pacific	+61 2 9938 8700
Other regions	+1 (425) 951-1330, or call your local representative
Fax	+1 (425) 951-6700
E-mail	General: ffss-service@fujifilm.com Europe, Middle East, and Africa: eraf-service@fujifilm.com Asia and Pacific: ffss-apacme-service@fujifilm.com
Web	www.sonosite.com


Getting started

All references to ECG have been removed from this chapter and the following sections shown.

System controls

The ECG control has been removed from the system controls.

- | | | |
|----|----------------------------|--|
| 1 | Control knobs | Turn to adjust gain, depth, cine buffer, brightness, and more, depending on context. Current functions appear on-screen above the knobs. |
| 2 | Freeze key | Press and hold to freeze or unfreeze the image. |
| 3 | Touchpad | When the touchpad is lit, use it to control items displayed on the screen. Double-tap the touchpad to switch between functions. |
| 4 | Touchpad key | Works in conjunction with the touchpad. Tap to activate an item on-screen, or to switch between functions. |
| 5 | Print key | Available only when a printer is connected to the system. Tap to print from a live or frozen scan. |
| 6 | Save keys | Tap one of these keys to save an image or a clip. |
| 7 | Image mode | Tap one of these keys to change the imaging mode. |
| 8 | System controls | Change system settings, switch transducers, add labels, or see patient information. |
| 9 | Image and Doppler controls | Use these to adjust the image or select the Doppler imaging mode. |
| 10 | Touchscreen | Use the touchscreen the same way you would use the touchpad. |

Intended uses

ECG has been removed from the following intended use.

Cardiac imaging applications

You can assess the heart size and function, cardiac valves, great vessels, visualize blood flow through cardiac valves, and assess for the presence or absence of pathology. In addition, you can identify the presence and location of fluid around the heart and lungs used to assist in pericardiocentesis and thoracentesis procedures. You can detect normal lung motion for the presence or absence of pathology.

Imaging

2D imaging

ECG has been removed from the 2D controls.

Table 4-1: 2D controls

Control	Description
Guide	Turns needle guidelines on. Guidelines can be used for needle guidance, and depend on transducer type. If using a variable angle needle guide, tap Guide . To select the angle, tap A , B , or C . To change the depth, move your finger on the touchscreen or the touchpad. To turn needle guidelines off, tap A , B , or C until the word Guide appears.

ECG

The ECG section has been removed from the Imaging chapter.

Measurements and calculations

Cardiac calculations

The ECG measurement has been removed as one of the ways to obtain an HR value from the list below. The section with details about the ECG heart rate measurement below the list has also been removed.

When performing cardiac calculations, the system uses the heart rate (HR) value present in the patient information form. The HR value can be obtained in any three different ways:

- ▶ Manual entry in the patient information form

- ▶ Doppler measurement
- ▶ M-Mode measurement

Cleaning and disinfecting

Cleaning and disinfecting the ECG cable and slave cable

This section has been removed from this chapter.

Safety

Electrical safety classification

Type CF applied parts With ECG module/ECG leads has been removed from the following list.

Class I equipment

The ultrasound system is classified as Class I equipment when powered from the external power supply or mounted on the stand because the external power supply is a Class 1 protectively earthed power supply. The stand has no protective earth. Ground bond testing is not applicable to the ultrasound system or the stand.

Note: AC powered peripherals that may be used with the system are Class I and are individually protectively earthed. Ground bond testing may be conducted on each AC powered peripheral.

Internally powered equipment

Ultrasound system not connected to the power supply (battery only)

Type BF applied parts

Ultrasound transducers

IPX-7 (watertight equipment)

Ultrasound transducers

IPX-8 (watertight equipment)

Footswitch

Non AP/APG

Ultrasound system power supply, SonoSite SII stand, and peripherals. Equipment is not suitable for use in the presence of flammable anaesthetics.

Electrical safety

ECG has been removed from the following warnings.

WARNINGS

- To avoid the risk of electrical shock:
- ▶ Connect this equipment to a supply main with protective earth.
 - ▶ Do not allow any part of the system (including the bar code scanner, external mouse, power supply, power supply connector, external keyboard, and so on), to touch the patient.
 - ▶ Use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade" or the equivalent. Do not remove or defeat the grounding wire.
 - ▶ When using the system in an environment where the integrity of the protective earth conductor arrangement is in doubt, operate the system on battery power only, without using the power supply.
 - ▶ Do not touch the power supply and the patient at the same time.
 - ▶ Do not touch any of the following:
 - ▶ The signal input/output connectors on the back of the ultrasound system
 - ▶ The system battery contacts (inside the battery compartment)
 - ▶ The system transducer connectors when the transducers are disconnected
 - ▶ Do not connect either of the following to an MSO or extension cord:
 - ▶ System power supply
 - ▶ Auxiliary mains outlet receptacles on the SonoSite SII stand

Compatible accessories and peripherals

ECG has been removed from the following table.

Table 9-2: Accessories and peripherals compatible with SonoSite SII ultrasound system

Description	Maximum Cable Length ^a
Black & white printer	—
Black & white printer power cable	3.3 ft/1 m
Black & white printer control cable	6.0 ft/1.8 m
Black & white printer video cable	6.2 ft/1.9 m
Footswitch	9.8 ft/3 m

Footswitch USB extension cable	6.5 ft/2 m
SonoSite SII Stand	—
Power cord (system)	10 ft/3 m
Power supply with DC cable	6.8 ft/2 m
Power supply AC cable	39 in/1 m
PowerPark	—
USB wireless adapter	—
USB memory drive	—

^aFor transducers, the maximum cable length is measured between the strain reliefs. The stated length do not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.

^bThe L52x transducer is for veterinary use only.

Labeling symbols

The China Pollution Control (5) symbol has been added to the Standard labeling symbols table 9-6.

Table 9-6: Standard labeling symbols

Symbol	Title	Standards development organization	Reference number	Description
	China Pollution Control (5)	Ministry of Industry and Information Technology	—	Pollution Control Logo. (Applies to all parts/products listed in the China RoHS disclosure table. May not appear on the exterior of some parts/products because of space limitations.)

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